

Airstream Max Class II, Microbiological Safety Cabinet, Model AC2-4L, with optional support stand.

Airstream[®]MAX[®]

Class II, Microbiological Safety Cabinets
The Industry's Most Comfortable, Well-Lit Cabinet



ESCO

WORLD CLASS. WORLDWIDE.




Main Features

- Unique Esco Dynamic Chamber™ plenum design delivers quiet, uniform airflow. Negative pressure plenum surrounds contaminated positive pressure plenum; no fabric bags are used.
- Dual, long-life ULPA filters (per IEST-RP-CC001.3), for supply and exhaust airflow. Sentinel™ Gold microprocessor monitors airflow in real-time via temperature compensated sensors and conserves UV lamp life with unique timer function.
- Ergonomically angled front improves reach and comfort.
- Frameless, shatterproof sash is easier to clean, offers larger, unobstructed viewing area.
- Multi-piece tray components which lift and remove to provide easy access encourage surface decontamination and are autoclavable.
- The front sash is motorized for convenient one touch operation.
- Raised airflow grille maintains safety by preventing blockage.
- Transparent side windows, angled front, and reduced noise levels combine to create the most comfortable, well-lit cabinet in Esco's range.
- The Accuflow microprocessor based speed controller maintains constant, stable airflow despite building supply voltage fluctuations.
- Esco **ISOCIDE™** antimicrobial surface on all painted surfaces minimizes contamination.

EN
12469

Airstream Max Class II Microbiological Safety Cabinet features glass sides to enhance visibility inside the work area. Model AC2-4L, with optional support stand.



Airstream•MAX

Biological Safety Cabinets • Class II Microbiological Safety Cabinets (L-Series)



Operator, Product and Environmental Protection

Esco Airstream Max Class II Microbiological Safety Cabinets offer a premium level of operator, product and environmental protection with advanced technology at an economical price. Airstream Max Microbiological Safety Cabinets provide protection against Biosafety Levels 1, 2 and 3, and can be used for handling Biosafety Level 4, provided that the operator wears a positive pressure suit.

Containment and Protection

- A combination of a supply ULPA filter and an exhaust ULPA filter creates a fully integrated performance envelope for product, operator and environmental protection.
- Inflow of room air enters the front air grille to establish operator protection; room air does not enter the work zone, preventing product contamination.
- Raised armrest prevents the likelihood of inflow grille blocking by operator's arms.
- The inflow velocity, downflow velocity, air flow path and intake geometry are precision tuned and tested to create an optimum air curtain on the front aperture. This curtain maintains operator and product protection even in the unlikely event of a severe inflow or downflow imbalance that would compromise protection in a conventional cabinet.

Integrated Filtration System

Independent supply and exhaust filters provide 99.999% typical efficiency for particle sizes of 0.1 to 0.3 microns. Airstream Max filters meet the IEST-RP-CC001.3 recommended practice for ULPA performance (USA), and EN 1822 for H14 performance (EU).

- ULPA filters (per IEST-RP-CC001.3), are tested to a typical efficiency of >99.999% for 0.1 to 0.3 micron particles; these provide better filtration capability than conventional H13 HEPA filters that have a typical efficiency of > 99.99% for 0.3 micron particles.
- Modern separatorless mini-pleat filter construction maximizes the filter surface area to extend filter life and eliminate possible filter media damage by thin and sharp aluminum separators used in conventional HEPA filter construction.

Mini-pleat Separatorless Filter (left) vs. Conventional Aluminum Separator Filter (right)



Esco cabinets use Swedish Camfil Farr® mini-pleat filters without aluminum separators to increase filter efficiency, minimize the chance of leakage, and to prolong filter life. Filters include a lightweight aluminum frame for structural stability and elimination of swelling common to conventional wood frames.

- The filter assembly is constructed in accordance with EN1822 requirements for performance and fire retardant properties.
- The supply filter provides ISO Class 3 (per ISO14644.1) clean air to the work surface in a gentle vertical laminar flow for product protection.
- The exhaust filter traps biohazard particles acquired from the work surface before air is exhausted to the room, offering operator and environmental protection.

- The exhaust filter media is protected from mechanical damage by an integrated metal screen guard, which is absent from conventional HEPA filters.

Blower System

The AC2-L blower system is designed for high performance operation, maximum energy efficiency and minimal maintenance.


- Improved energy efficiency lowers operating costs.
- Reduced noise and vibration levels compared to conventional blowers provide a comfortable working environment.
- An integral blower hour meter tracks operating life and aids in predictive maintenance planning.
- The external rotor motor design allows for optimum cooling of the motor during extended operations and extends the motor bearing life.
- To prevent fan damage, a paper-catch grille traps papers or towels that may drop down on the drain pan, preventing them from being pulled into the column by fan suction.

Sentinel™ Microprocessor Control, Alarm, Monitoring System


The Esco Sentinel™ microprocessor-based control system supervises operation of all cabinet functions.

- The control panel is located on the center of the cabinet and is angled down for easy access by the operator.
- Continuous monitoring of cabinet airflow is displayed on a bright, easy-to-read LCD panel. The large display monitors operational parameters.
- True temperature-compensated airflow velocity sensors provide independent

Esco Accuflow Microprocessor Speed Controller



The Esco Accuflow™ microprocessor speed controller maintains steady motor/blower speed despite building voltage fluctuations, thereby assuring constant face velocity and downflow for optimum safety, containment and protection.



Cabinet Filtration System

Dynamic air barrier, inflow and forward-directed downflow air converge

- Ambient air is pulled through the perforations located towards the work zone front to prevent contamination of the work surface and work product. The inflow does not mix with the clean air within the cabinet work zone. Inflow air travels through a return path toward the common air plenum (blower plenum) at the top of the cabinet.
- The uniform, non-turbulent air stream protects against cross contamination within and throughout the work area.
- Near the work surface, the downflow air stream splits with a portion moving toward the front air grille, and the

remainder moving to the rear air grille. A small portion of the ULPA filtered downflow enters the intake perforations at the side capture zones at a higher velocity (small blue arrows).

- A combination of inflow and downflow air streams forms an air barrier that prevents contaminated room air from entering the work zone, and prevents work surface emissions from escaping the work zone.
- Air returns to the common air plenum where the 33% exhaust and 67% recirculation process is continued.

- ULPA-filtered air
- Unfiltered / potentially contaminated air
- Room air / Inflow air

4

measurement of inflow and downflow velocities despite room temperature fluctuation.

- All electronic parts are contained inside a plug-and-play module that permits easy exchange if required.
- Microprocessor software updates are available from Esco for download via the Internet.

Consult your Esco Operating Manual or contact Esco or your Sales Representative for information on user-preference programming capabilities built into the Sentinel microprocessor platform.

Cabinet Construction

Robust construction and enhanced safety features qualify the cabinet for the most demanding laboratory applications. The cabinet is fully assembled and ready to install and operate when shipped.

- The cabinet work zone has no welded joints to collect contaminants or rust.
- All stainless steel work surfaces are accessible for cleaning.

- Work zone back wall is manufactured from powder coated steel, which eliminates glare associated with stainless steel.
- Multi-piece tray components lift and remove to provide easy access and to encourage surface decontamination.
- A recessed central area and drain pan channels spills and prevents liquids from entering the lower filtration and blower systems.
- There are no screws on the front or sides to trap contaminants or complicate cleaning.
- External and internal surfaces are coated with Esco Isocide™ antimicrobial coating to protect against surface contamination and inhibit bacterial growth. Isocide eliminates 99.9% of surface bacteria within 24 hours of exposure.

Service Fitting Access

The cabinet is prepared for easy installation of optional gas and vacuum fittings; see Accessories.

- Optional service fittings openings are offset for easier access.

Comfortable Ergonomic Design


The cabinet is engineered for comfort, utility value and safety.

- The angled viewing window and narrow profile front grille improve reach into the work area.
- The instant-start 5000k fluorescent lamp operates on an electronic ballast to reduce heat, improve comfort and conserve energy.
- The lamp delivers uniform lighting to the work surface for greater comfort, reduced glare and improved productivity; see Specifications.
- The front armrest is raised above the work zone to improve comfort and to minimize blockage of forward airflow perforations.
- The optional adjustable support stand provides work surface height control.
- The frameless sash eliminates operator's line of sight blockage.
- A generous sash opening allows for easier access into the work zone, provides ample room for transferring of small equipment; see Specifications.
- The motorized window can be fully opened to insert and remove larger instrumentation and equipment.

Electrical Safety and Certification

All components meet or exceed applicable safety requirements.

Esco ULPA Filter Efficiency




Independent supply and exhaust filters provide 99.999% typical efficiency for particle sizes of 0.1 to 0.3 microns. Airstream Series filters meet the IEST-RP-CC001.3 recommended practice for ULPA performance (USA), and EN 1822 for H14 performance (EU).

Enlarged, multi-line digital read-out with alpha-numeric display indicates all input, status and alarm functions.

Enlarged touchpad data entry buttons with tactile feedback permit control settings and access to diagnostics, default settings and hierarchical menus.

Color coded indicator lamps display green for fan operation; blue for fluorescent lights and electrical outlet; and orange for UV lamp ON caution.

Programmable automatic UV light timer simplifies operation, enhances contamination control, extends UV lamp life and saves energy.



Sentinel Microprocessor Control System, Programmable

- When programmed ON
 - the start-up sequence confirms status with Air Safe and local time display.
 - the Personal Identification Number (PIN) access restricts unauthorized adjustments.
 - an airflow alarm warns of deviations from normal velocities.

- Each cabinet is individually factory tested for electrical safety.
- Documentation specific to each cabinet serial number is maintained on file.
- Tested to major world standards for microbiological safety cabinets, including the EN 12469.
- Contact EscO or your Sales Representative for site preparation information; see Electrical Specifications.

Warranty

Airstream Max AC2-L Series cabinets are warranted for 3 years excluding consumable parts and accessories.

- Each cabinet is shipped with a comprehensive User's Manual complete with a report documenting all test procedures.

- Additional IQ/OQ documentation is available upon request.
- Contact your local Sales Representative for specific warranty details or documentation requests.

Accessories and Options

EscO offers a variety of options and accessories to meet local applications. Contact EscO or your local Sales Representative for ordering information.

Support Stands


- Fixed height with casters (SPC)
 - Available 711 mm (28") or 860 mm (34")
- Fixed height with leveling feet (SAL)
 - Available 737 mm (29") or 838 mm (33"), ±38.1 mm (1.5")

- Hydraulic motorized adjustable height with casters (SPM)
 - Nominal range 711 mm (28") to 864 mm (34")
- Telescoping height with casters (STC)
 - Nominal range 660 mm (26") to 880 mm (34.6")
 - Adjustable in 25.4 mm (1") increments
- Telescoping height with leveling feet (STL)
 - Nominal range 660 mm (26") to 960 mm (37.8")

Electrical Outlets and Utility Fixtures

- Electrical outlet, ground fault, North America
- Electrical outlet, Europe / Worldwide

Robust Cabinet Construction and Enhanced Safety Features



Service fixtures are offset for easier reach. Standard cabinets include two fixture provisions on each sidewall (one provision on each sidewall for 0.9 meter/3 ft. cabinet).

- Helpful for certifiers, the hinged maintenance assembly opens to a fixed position on integrated, stainless steel struts providing front service access.



All key components, with the exception of the blower/motor assembly, are mounted outside the air stream and away from contaminated air to permit service without decontamination. These include fluorescent lamps, electrical harnesses, electronic boards and microprocessor control.

- Panels enclosing microbiologically or electrically hazardous areas are color-coded red to warn service technicians.
- The telescoping Dynamic Chamber™ plenum minimizes physical lifting and accelerates filter change when required.
- Work area containment is maintained even when removable components are lifted out for cleaning.



The multi-piece stainless steel work tray edges are radiused and easy to clean without crevices or joints.

- The lower drain trough is a single-piece fabrication with wide open angles and a channel to direct spills to the drain.
- The closed sidewall contains no perforations, air return slots or other hidden areas where contaminants can accumulate.

Dynamic Chamber™ Plenum Design



■ Negative pressure ■ Positive pressure


The Esco double-wall design creates a Dynamic Chamber plenum which surrounds contaminated areas with negative pressure, preventing the possibility of contamination from leaks in filter seal, gasket or cabinet structure.

- Petcock (air, gas, vacuum)
 - North America (American) style
 - Europe / Worldwide style DIN 12898, DIN 12919, DIN 3537

Cabinet Accessories

- Germicidal UV lamp
 - Controlled by automatic UV lamp timer through Sentinel™ micro-processor control panel
 - Emission of 253.7 nanometers for most efficient decontamination
- PVC armrest
 - Chemically treated, improves operator comfort, easy to clean.
- Ergonomic lab chair
 - Laboratory grade construction, meets Class 100 cleanliness; alcohol resistant PVC materials
 - Adjustable height 395-490 mm (15.6" - 19.3")
- Ergonomic foot rest
 - Angled, helps maintain proper posture
 - Adjustable height
 - Anti-skid coating, chemical resistant finish
- IV bar, with hooks
 - Stainless steel construction
 - Available for all standard Esco cabinets
- Microscope viewing device
 - Mounting and viewing pouch integrated into sash. Factory installed; specify when ordering.

Airstream Max Model AC2 (L-Series), Class II Microbiological Safety Cabinet Technical Specifications



1. Exhaust ULPA filter
2. Blower
3. Downflow ULPA filter
4. Downflow sensor
5. Standard UV light Retrofit Kit provision

6. Plugged service fixtures provisions (2 on each side)
7. Electrical/ Electronics panel
8. Fluorescent lamps
9. Tempered side glass
10. Stainless steel arm rest

11. Exhaust sensor
12. Esco Sentinel microprocessor control system
13. Motorized sash window
14. Female IEC inlet for UV provision
15. Standard electrical outlet Retrofit Kit provision
16. Stainless steel multi-piece work tray

General Specifications, Airstream Max Class II, Microbiological Safety Cabinets (L-Series)

Model		AC2-4L1	AC2-6L1
Nominal Size		1.2 meters (4')	1.8 meters (6')
External Dimensions (W x D x H)	Without Base Stand	1340 x 740 x 1413 mm 52.8" x 29.1" x 55.6"	1950 x 740 x 1413 mm 76.7" x 29.1" x 55.6"
	With Optional Base Stand, 711 mm (28") type	1340 x 740 x 2124 mm 52.8" x 29.1" x 83.6"	1950 x 740 x 2124 mm 76.7" x 29.1" x 83.6"
Internal Work Area, Dimensions (W x D x H)		1270 x 560 x 670 mm 50.0" x 22.0" x 26.4"	1880 x 560 x 670 mm 74.0" x 22.0" x 26.4"
Internal Work Area		0.58 m ² (6.2 sq.ft)	0.87 m ² (9.3 sq.ft)
Tested Opening		200 mm (7.9")	200 mm (7.9")
Working Opening		225 mm (8.9")	225 mm (8.9")
Average Airflow Velocity	Inflow	0.45 m/s (90 fpm) at initial setpoint	
	Downflow	0.30 m/s (60 fpm) at initial setpoint with uniformity of better than +/- 20%	
Airflow Volume	Inflow	412 m ³ /h (245 cfm)	606 m ³ /h (356 cfm)
	Downflow	741 m ³ /h (436 cfm)	1096 m ³ /h (645 cfm)
	Exhaust	412 m ³ /h (245 cfm)	606 m ³ /h (356 cfm)
ULPA Filter Typical Efficiency	Downflow	>99.999% at 0.1 to 0.3 microns and MPPS as per IEST-RP-CC001.3 USA with H14 rating as per EN 1822, Europe	
	Exhaust		
Typical Sound Emission per EN 12469*		<55 dBA	<57 dBA
Fluorescent Light Intensity At Zero Ambient		>1250 Lux (>116 foot candles)	1250 Lux (>116 foot candles)
Cabinet Construction	Main Body	1.2 mm (0.05") 18 gauge electrogalvanized steel with white oven-baked epoxy-polyester Isocide antimicrobial powder coated finish	
	Work Surface	1.5 mm (0.06") 16 gauge stainless steel, type 304, with BA finish	
	Side Walls	UV absorbing tempered glass, 5 mm (0.2"), colorless and transparent	
Electrical 220-240V, AC, 50Hz, 1Ø	Cabinet Power/ Amp	550 W/ 2.5 A	650 W/ 3.5 A
	Outlet Amp Fuse	5 A	5 A
	Full Load Amps	7.5 A	8.5 A
	BTU/ Hr	1122	1326
Net Weight**		258 kg (585 lbs)	300 kg (661 lbs)
Shipping Weight**		313 kg (710 lbs)	354 kg (780 lbs)
Shipping Dimensions, Maximum (W x D x H)**		1500 x 950 x 1800 mm 59.0" x 37.4" x 70.8"	2100 x 950 x 1800 mm 82.7" x 37.4" x 70.8"
Shipping Volume, Maximum**		2.57 m ³ (91 cu.ft.)	3.59 m ³ (127 cu.ft.)

* Noise reading in open field condition/ anechoic chamber.

** Cabinet only; excludes optional stand

Standards Compliance	For Microbiological Safety Cabinets	For Air Quality	For Filtration	For Electrical Safety
	EN 12469, Europe	ISO 14644.1 Class 3, Worldwide AS 1386 Class 1.5, Australia JIS B9920 Class 3, Japan	EN-1822 (H14), Europe IEST-RP-CC001.3, Worldwide IEST-RP-CC007.1, Worldwide IEST-RP-CC034.1, Worldwide	IEC 61010-1, Worldwide EN 61010-1, Europe UL 61010-1, USA CAN/CSA-22.2, No.61010-1



- Biological Safety Products
- Cleanroom Products
- Containment / Pharma Products
- Ductless Fume Hoods / Carbon Filtration
- General Purpose Scientific Equipment
- Industrial Lab Equipment
- In-Vitro Fertilization Products
- PCR Products
- Pharmacy Products
- Lab Animal Research Products
- Lab Thermostatics Products
- Lab Ventilation / Chemical Fume Products / Lab Furniture
- Powder Handling Products

Since 1978, Esco has emerged as a leader in the development of controlled environment, laboratory and cleanroom equipment solutions. Products sold in more than 100 countries include biological safety cabinets, fume hoods, ductless fume hoods, laminar flow clean benches, animal containment workstations, cytotoxic cabinets, hospital pharmacy isolators, and PCR cabinets and instrumentation. With the most extensive product line in the industry, Esco has passed more tests, in more languages, for more certifications, throughout more countries than any biosafety cabinet manufacturer in the world. Esco remains dedicated to delivering innovative solutions for the clinical, life science, research and industrial laboratory community. www.escoglobal.com.

Biological Safety Cabinets and Laminar Flow • Laboratory Fume Hoods • Laboratory Ovens
 Laboratory Incubators • PCR Thermal Cyclers • Microplate Shaker/Incubators • Ultraflow Freezers

ESCO

WORLD CLASS. WORLDWIDE.

Esco Technologies, Inc. • 2940 Turnpike Drive, Units 15-16 • Hatboro, PA 19040, USA
 Toll-Free USA and Canada 877-479-ESCO • Tel 215-441-9661 • Fax 215-441-9660
us.escoglobal.com • usa@escoglobal.com

Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777
 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com
www.escoglobal.com

Esco Global Offices | Kuala Lumpur, Malaysia | Leiden, The Netherlands | Manama, Bahrain
 Mumbai, India | Philadelphia, USA | Salisbury, UK | Shanghai, China | Singapore

