



ColdFAST

DW-HL778M



Technical specifications

ColdFAST

DW-HL778M

DW-HL778M

| | | |
|-----------------|-----------------------------------|---|
| Model | Type of Construction | Upright |
| | Capacity (l) | 778 |
| | Inner Dimensions (W x D x H) mm | 865 x 696 x 1286 |
| | Outer Dimensions (W x D x H) mm | 1205 x 1025 x 1955 |
| | Package Dimensions (W x D x H) mm | 1320 x 1155 x 2171 |
| | Net / Gross Weight (kgs) | 359 / 403 |
| Performance | Cooling Performance (°C) | -86 |
| | Temperature Range (°C) | -40 ~ -86 |
| | Ambient Temperature (°C) | 16 ~ 32 |
| | Climate Class | N |
| | Sound Level (dB(A)) | 51 |
| | Controller | Microprocessor |
| Refrigeration | Display | 10" HD Intelligent Touch Screen |
| | Compressor | 2 pcs |
| | Cooling Mode | Direct Cooling / Single Cascade |
| | Defrost Mode | Manual |
| | Refrigerant | HC |
| | Insulation Thickness (mm) | 130 |
| Construction | External Material | Cold Rolled Steel Sheet |
| | Inner Material | Stainless Steel |
| | Shelves | 3 (Stainless Steel) |
| | Inner doors | 2 (optional: 4) |
| | Door Lock | Yes, with key |
| | Padlock | Yes |
| | Access Test Port | 3 pcs. Ø 25 mm |
| | Casters | 4 + (2 leveling feet) |
| | Data Logging / Time / Duration | USB / Record Every 5 Minutes / 10 Years |
| | Backup Battery | Yes |
| Alarm | WiFi | Yes |
| | Temperature | High / Low Temperature, High Ambient Temperature |
| | Electrical | Power Failure, Low Battery |
| Electrical Data | System | Sensor Failure, Main Board Communication Error, Condenser Overheating, Door Ajar |
| | Power Supply (V/Hz) | 220-240~ / 50 |
| | Rated Power (W) | 1380 |
| | Current Consumption (KWh/24h) | 10,0 |
| Accessories | Electrical Current (A) | 11,6 |
| | Standard | RS485, Remote Alarm Contact |
| | Options | Chart-Recorder, CO ₂ -Backup-System, Thermal-Printer, Bluetooth-Printer, Finger-Print or RFID-Chip, Racks, Boxes & Tubes |

The model, parameters and performance specified in this brochure may be changed without prior notice because of product upgrading. There may be differences between the product images shown in this brochure and the actual products. When you are buying any product, please check the actual product.