



# ERBAdry<sup>®</sup>: ANHYDROUS SOLVENTS

## *new generation*

Our innovation lies in the caps of our anhydrous ERBAdry<sup>®</sup> solvents. Indeed, these new generation caps and septa have been developed by CARLO ERBA Reagents for improved ergonomics and better water content retention over time:

- 3 large areas available to improve usability
- Septum double layers to prevent chemical from contamination
- 2 types of septa selected for better chemical compatibility
- Flip Top Cap always ready to be closed
- 2 temper evident rings showing the double unused product warranty
- 100ml & 1l amber glass bottles to protect chemicals from UV



**AN INNOVATION  
CARLO ERBA Reagents**



**CARLO ERBA**  
REAGENTS



# ERBAdry®: ANHYDROUS SOLVENTS

*new generation*

## AN INNOVATION CARLO ERBA Reagents

Our first selection of solvents available from November (table below). Others, planned very soon, will be announced to you on our website or through our Newsletters.

Product	Pkg	Code
Acetone RS ERBAdry®	100 ml	P00510D03
	1 l	P00510D16
Acetonitrile RS ERBAdry®	100 ml	P00610D03
	1 l	P00610D16
Chloroforme stabilized with amylene RS ERBAdry®	100 ml	P02410AD03
	1 l	P02410AD16
Chloroforme stabilized with ethanol RS ERBAdry®	100 ml	P02410ED03
	1 l	P02410ED16
Cyclohexane RS ERBAdry®	100 ml	P02510D03
	1 l	P02510D16
Dichloromethane stabilized with amylene RS ERBAdry®	100 ml	P02910AD03
	1 l	P02910AD16
Diethyl ether RS ERBAdry®	100 ml	P013A10D03
	1 l	P013A10D16
n,n-Dimethylformamide RS ERBAdry®	100 ml	P03410D03
	1 l	P03410D16
1,4-Dioxane RS ERBAdry®	100 ml	P03610D03
	1 l	P03610D16
Ethanol absolu RS ERBAdry®	100 ml	P013A10D03
	1 l	P013A10D16
Ethyle acetate RS ERBAdry®	100 ml	P00210D03
	1 l	P00210D16
Isopropyl ether RS ERBAdry®	100 ml	P04310D03
	1 l	P04310D16
Methanol RS ERBAdry®	100 ml	P09310D03
	1 l	P09310D16
Propanol-2 RS ERBAdry®	100 ml	P09510D03
	1 l	P09510D16
Pyridine RS ERBAdry®	100 ml	P06710D03
	1 l	P06710D16
Toluene RS ERBAdry®	100 ml	P07110D03
	1 l	P07110D16



Orange septum : Butyl rubber / PTFE



Black septum : Bromobutyl / PTFE

CER FR/ ERBADRY® : SOLVANTS ANHYDRES/2021-10/Ed 01. All pictures and specifications included in this document are purely indicative and may be subject to change without notice.



WWW.CARLOERBAREAGENTS.COM

"Partner in your choice"

CARLO ERBA Reagents S.A.S.  
CHAUSSEE DU VEXIN  
27106 VAL DE REUIL  
TEL. +33 2 32 09 20 00  
FAX +33 2 32 09 20 20  
export@cer.dgroup.it  
www.carloerbareagents.com

