



SmartFAST Light

Technical data sheet



Technical specifications

FASTER product protection workbench

FASTER SmartFAST Mini | ISO Class 3 according to DIN EN ISO 14644-1

Description	Unit	0.6 m
General dimensions W×H×D	mm	642×830×330
Usable dimensions W×H×D	mm	580×375×317
Working height (seated workstation)	mm	865
Work opening	mm	200
Weight	kg	50
Laminar flow	m/s	0,4
Noise level	dB(A)	<60
Electrical connection		1Ph+E - 230 V 50 Hz
Power consumption	W	130
Protection class / IP		1/20
Internal fuse protection	kW	0,11

CER DE/Tec_SmartFAST-Mini_EN/2022-08/Ed. 01 In this brochure all pictures and specifications are purely indicative and may vary without notice.