



SmartFAST

Technical data sheet



Technical specifications

FASTER product protection workbench

FASTER SmartFAST | ISO Class 3 according to DIN EN ISO 14644-1

Description	Unit	0.8 m
General dimensions W×H×D	mm	830×910×652
Usable dimensions W×H×D	mm	730×600×378
Unit height with underframe (standard)	mm	1720
Working height (seated workstation)	mm	800
Work opening	mm	200
Weight	kg	85
Laminar flow	m/s	0,4
Noise level	dB(A)	<60
Light intensity	lux	565
Electrical connection		1Ph+E - 230 V 50 Hz
Power consumption	W	178
Protection class / IP		1/20
Internal fuse protection	A	4
Sockets		1
Other sockets		●
Aluminium front blind		●
Automatic flow regulator		●

● Standard ● Optional ○ Not available

CER DE/Tec_SmartFAST_EN/2022-08/Ecl.01 In this brochure all pictures and specifications are purely indicative and may vary without notice.